

Chemistry 101
MOLES, CHEMICAL FORMULAS, AND STOICHIOMETRY
Pre-Lab Exercises

Student: _____

Date: _____

Instructor: _____

Section: _____

1. Calculate the molar mass of $\text{CoCl}_2 \cdot 6\text{H}_2\text{O}$.

g/mol

2. Calculate the percent by mass of H_2O in $\text{CoCl}_2 \cdot 6\text{H}_2\text{O}$.

%

3. Balance the equation:



4. How many moles are in 4.56 g of KClO_3 ?

mol

5. Based upon the balanced chemical equation from #3, how many grams of KCl will form when 4.56 g of KClO_3 decompose?

g

Chemistry 101


MOLES, CHEMICAL FORMULAS, AND STOICHIOMETRY

Student: _____
Partner: _____
Instructor: _____
Section: _____ Date: _____

PROCEDURE AND REPORT

Formulas and Dehydration of Hydrates

Hydrates form when water molecules are bound within the crystal structure of a solid compound. A common example of a hydrate is $\text{CuSO}_4 \cdot 5\text{H}_2\text{O}$. In this case, 5 moles of water are bound for every one mole of CuSO_4 in the crystal structure. The name of this hydrate is copper (II) sulfate pentahydrate. The water molecules are written separated from the rest of the formula by a dot because they are bound to the CuSO_4 compound as regular, distinct water molecules. Although hydrates are quite stable, the water molecules typically can be removed from the crystal structure by applying heat. This process is called **dehydration** and the compound remaining is called an **anhydrous salt**. In the case of $\text{CuSO}_4 \cdot 5\text{H}_2\text{O}$, heating the hydrate to a high enough temperature leads to the anhydrous salt CuSO_4 and releases H_2O as indicated by the equation:



When calculating the **molar mass** of a hydrate, the mass of the water molecules must be included. The molar mass of CuSO_4 is 159.61 g/mol. The molar mass of $\text{CuSO}_4 \cdot 5\text{H}_2\text{O}$ is 249.69 g/mol. In one mole of $\text{CuSO}_4 \cdot 5\text{H}_2\text{O}$, there are 159.61 g CuSO_4 and 90.1 g (5×18.02 g) of H_2O . To find the percent by mass of water in $\text{CuSO}_4 \cdot 5\text{H}_2\text{O}$, simply take the ratio of the mass of H_2O in $\text{CuSO}_4 \cdot 5\text{H}_2\text{O}$ to the mass of $\text{CuSO}_4 \cdot 5\text{H}_2\text{O}$ times 100%:

$$\% \text{water} = \frac{\text{mass } \text{H}_2\text{O}}{\text{mass } \text{CuSO}_4 \cdot 5\text{H}_2\text{O}} \times 100\% = \frac{90.1\text{g}}{249.69\text{g}} \times 100\% = 36.1\%$$

Dehydration of a Known Hydrate


Start the experimental procedure outlined below. However, during periods of time when you are waiting for your crucible to cool, perform the following preliminary calculations.

Preliminary Calculations

- Calculate the molar mass of MgSO_4 : _____
- Calculate the molar mass of H_2O : _____
- Calculate the molar mass of $\text{MgSO}_4 \cdot 7\text{H}_2\text{O}$: _____
- Based upon the formula, calculate the percent, by mass, of MgSO_4 in $\text{MgSO}_4 \cdot 7\text{H}_2\text{O}$: _____
- Based upon the formula, calculate the percent, by mass, of H_2O in $\text{MgSO}_4 \cdot 7\text{H}_2\text{O}$: _____

Dehydration of a Known Hydrate (continued)

- Clean and dry a small crucible. Support the crucible using a ringstand and wire triangle. Heat the crucible over a Bunsen burner, gently for a few seconds then strongly for 2 minutes. Turn off the flame and allow 5 minutes for the crucible to cool. Once the crucible has been heated, handle it with tongs to prevent moisture and oil from your hands from changing the mass of the crucible. As long as you touch the crucible only with tongs, you will not have to repeat this step in future trials.
- When the crucible has cooled, measure and record the mass of the crucible to the nearest 0.001 g. **Never place a hot object on the balance!!**
- Fill the crucible about $\frac{1}{4}$ full of $\text{MgSO}_4 \cdot 7\text{H}_2\text{O}$. Measure the mass of the crucible and $\text{MgSO}_4 \cdot 7\text{H}_2\text{O}$ to the nearest 0.001 g.
- Calculate the mass of $\text{MgSO}_4 \cdot 7\text{H}_2\text{O}$ in the crucible by subtraction.
- Dehydrate the $\text{MgSO}_4 \cdot 7\text{H}_2\text{O}$ by heating the crucible over a Bunsen burner, gently for about 1 minute then strongly for 5 minutes. Allow the crucible to cool for 5 minutes.
- When the crucible has cooled, measure and record the mass of the crucible and the anhydrous MgSO_4 to the nearest 0.001 g. Discard the crucible contents. If you do not touch the crucible with your hands, you will not need to clean the crucible by heating before future trials.
- Calculate the mass of dehydrated MgSO_4 in the crucible after heating.
- Calculate the mass of H_2O released from the hydrate during heating.
- Calculate the percent by mass of H_2O in the original hydrate (mass of H_2O /Mass of $\text{MgSO}_4 \cdot 7\text{H}_2\text{O}$ x 100%). How closely do your experimental results correlate with the theoretical value calculated above in preliminary calculation e?
- Calculate the moles of dehydrated MgSO_4 in the crucible.
- Calculate the moles of H_2O released from the hydrate.
- Calculate the ratio of moles of H_2O to moles of MgSO_4 . Compare your result with the number you expect based upon the formula of $\text{MgSO}_4 \cdot 7\text{H}_2\text{O}$.
- Perform a second trial if instructed to do so by you instructor.

**Dehydration of a Known Hydrate**

Steps	Trial #1	Trial #2
2. Mass of crucible		
3. Mass of crucible and $\text{MgSO}_4 \cdot 7\text{H}_2\text{O}$		
4. mass of $\text{MgSO}_4 \cdot 7\text{H}_2\text{O}$ in crucible		
6. Mass of crucible and MgSO_4 after heating		
7. Mass of MgSO_4 remaining after heating		
8. Mass of H_2O released		
9. % H_2O in $\text{MgSO}_4 \cdot 7\text{H}_2\text{O}$		
10. Moles MgSO_4		
11. Moles H_2O released		
12a. Calculated Moles H_2O /Moles MgSO_4		
12b. Expected Moles H_2O /Moles MgSO_4		

Calculations:

Formula of an Unknown Hydrate

In this section you will determine the value of X in the formula of a hydrate of calcium sulfate, $\text{CaSO}_4 \cdot \text{XH}_2\text{O}$.

1. If necessary, clean a crucible as previously instructed.
2. When the crucible has cooled, measure and record the mass of the crucible to the nearest 0.001 g. **Never place a hot object on the balance!!**
3. Fill the crucible about $\frac{1}{4}$ full of $\text{CaSO}_4 \cdot \text{XH}_2\text{O}$. Measure the mass of the crucible and $\text{CaSO}_4 \cdot \text{XH}_2\text{O}$ to the nearest 0.001 g.
4. Calculate the mass of $\text{CaSO}_4 \cdot \text{XH}_2\text{O}$ in the crucible.
5. Dehydrate the $\text{CaSO}_4 \cdot \text{XH}_2\text{O}$ by heating the crucible over a Bunsen burner, gently for about 1 minute then strongly for 5 minutes. Allow the crucible to cool for 5 minutes.
6. When the crucible has cooled, measure and record the mass of the crucible and the anhydrous CaSO_4 to the nearest 0.001 g. Discard the crucible contents. If you do not touch the crucible with your hands, you will not need to clean the crucible by heating before future trials.
7. Calculate the mass of dehydrated CaSO_4 in the crucible after heating.
8. Calculate the mass of H_2O released from the hydrate during heating.
9. Calculate the percent by mass of H_2O in the original hydrate.
10. Calculate the moles of dehydrated CaSO_4 in the crucible.
11. Calculate the moles of H_2O released from the hydrate.
12. Find X in $\text{CaSO}_4 \cdot \text{XH}_2\text{O}$ by calculating the ratio of moles of H_2O to moles of CaSO_4 .
13. Round the value determined in step 12 to the nearest whole number and write the formula for the unknown hydrate.
14. Perform a second trial if instructed to do so by you instructor.

Formula of an Unknown Hydrate

Steps	Trial #1	Trial #2
2. Mass of crucible		
3. Mass of crucible and $\text{CaSO}_4 \cdot \text{XH}_2\text{O}$		
4. Mass of $\text{CaSO}_4 \cdot \text{XH}_2\text{O}$ in crucible		
6. Mass of crucible and CaSO_4 after heating		
7. Mass of CaSO_4 remaining after heating		
8. Mass of H_2O released		
9. % H_2O in $\text{CaSO}_4 \cdot \text{XH}_2\text{O}$		
10. Moles CaSO_4		
11. Moles H_2O released		
12. Moles H_2O /Moles CaSO_4		
13. Formula of unknown		

Calculations:



Decomposition of Baking Soda

When baking soda or sodium bicarbonate, NaHCO_3 , is heated, it decomposes releasing H_2O and CO_2 . The third decomposition product is a solid, either sodium carbonate, Na_2CO_3 , or sodium oxide, Na_2O . Balanced chemical equations for both of these reactions have been written below even though only one of the reactions occurs. In this section you will determine which of these two decomposition reactions actually occurs.

1. Use the crucible from the previous section. As long as you have not handled the crucible with your hands, it is not necessary to heat the crucible to clean it. Measure the mass of the crucible to the nearest 0.001g.
2. Add about 0.5 gram of baking soda to the crucible and record the new mass.
3. Calculate the mass of the baking soda in the crucible by subtraction.
4. Heat the crucible over a Bunsen burner, gently for about 1 minute then strongly for 2 minutes. Allow the crucible to cool for 5 minutes.
5. When the crucible has cooled, measure and record the mass of the crucible and decomposition product to the nearest 0.001 g. Discard the crucible contents.
6. Calculate the mass of the decomposition product.

1. Mass of crucible	
2. Mass of crucible and NaHCO_3	
3. Mass of NaHCO_3 in crucible	
5. Mass of crucible and product after heating	
6. Mass of product remaining after heating	

The chemical equations for the two possible decomposition reactions are shown below. Based upon the initial mass of NaHCO_3 , complete the following two tables to determine the expected mass of each of the potential products. You will then use the expected mass of each product and the measured product mass (line #6) to determine which reaction actually occurred.



Initial mass of NaHCO_3 (from line #3)	
Expected mass Na_2CO_3	

Initial mass of NaHCO_3 (from line #3)	
Expected mass Na_2O	

Compare the expected mass calculated for Na_2CO_3 and the expected mass calculated for Na_2O with the experimental mass of the product from line #6.

Write the formula for the decomposition product: _____

Write the correct decomposition reaction:

Calculations: